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©2011-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Stellant, MEDRAD Stellant, FluiDots, and MEDRAD FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

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ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for Use: The contents of this package are intended to be used to transfer contrast media into a syringe.

Contraindications: These devices are not intended for multiple patient use nor for any other use for which they are not indicated. Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

Air embolization can cause death or serious injury to the patient. Do not connect the patient to the injector until all trapped air has been cleared from the syringe and fluid path.

Biological contamination can result from reusing disposable items or failure to follow aseptic technique. Properly discard disposable items after use.

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used.

Patient or operator injury could result from using a closed stopcock during an injection. A closed stopcock could cause an overpressure condition causing the syringe to burst.

Patient or operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open: occlusions in the fluid path may result in leaks or ruptures to the transfer set.

Operator injury could result from using a sharp spike. Use care in handling and inserting the spike into the bottle.

The PVC tubing in this product contains DEHP, a chemical known to the State of California to cause birth defects or other reproductive harm.

Extremely high doses of phthalates such as DEHP have been reported to cause birth defects in research animals. Although limited scientific data is available on the effect of phthalates in humans, health authorities recommend that exposure to devices containing these compounds be minimized as much as possible for children and nursing or pregnant women.

DEHP is also a substance on the candidate list of Substances of Very High Concern (SVHC) according to Article 59 of REGULATION (EC) no. 1907/2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH).

Substance name: Bis (2-ethylhexyl)phthalate (DEHP)
CAS Number: 117-81-7
Ec Number: 204-211-0

Special Precautions: CVTS & MRTS

Patient or operator injury may result from damage to the dual check-valve. Do not use with flow rates higher than 4 milliliters per second.

Patient or operator injury may result from damage to the dual-check valve. Do not use with Angiographic Injectors or high pressure/dynamic applications greater than 300 PSI (20.7 BAR).

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten.

Refer to the Injector Operation Manual and syringe Instructions for Use for further instructions.

Transfer Set Installation

Table with 2 columns: Syringe Type, Number of Turns. Rows include FLS with Swivel Nut (1/2) and FLS with Threaded Luer (3/4).

Before installing the transfer set, ensure the syringe is attached to the injector and the injector head is tilted up.

CVTS and MRTS Installation

- 1. Connect the transfer set by inserting the contrast bottle into the holder, remove the spike cover and insert the spike into the bottle.
2. Remove the protective caps from the syringe and female luer of the transfer set and attach to the syringe tip.
3. Fill the syringe by retracting the syringe piston/plunger.

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